

Medical air Sterilizer & Purifier KJG200

Fitur :

- Desain yang unik
- HEPA Filter Mesh
- Active Carbon Filter Mesh
- Accelerant Deodorization Filter Mesh
- Panel Kontrol dengan 3 level pengaturan aliran udara



CPAKB
CARA PEMBUATAN
ALAT KESEHATAN
YANG BAIK
HK.02.04/VI/011/2016
PT. SINKO PRIMA ALLOY

Medical air Sterilizer & Purifier KJG200

SPESIFIKASI

Model		KJG200
Power Source		220-240V
Frequency		50/60 Hz
Output Power (W)		45W
Rated Flow Volume (m3/h)		150 m3/h
Application Rate (m2)		40 m2
Net Weight (kg)		± 5.2 KG
Outside Sizes (LxWxH mm)		370x180x500 mm
Kondisi Pemakaian Normal	Suhu	18°C ~ 30°C
	Kelembapan Relatif	10% ~ 85%
	Tekanan Atmosfir	700hpa ~ 1060hpa
Transportasi dan Penyimpanan	Suhu	-10°C ~ +50°C
	Kelembapan Relatif	≤ 95%
	Tekanan Atmosfir	500hpa ~ 1060hpa



84778/A/0001/UK/En



ISO 9001
REGISTERED FIRM

ISO9001:2015



globalgroup
ISO 13485

ISO13485:2016
MS 181050 R